Blood Pressure Monitor

For more information see the extended manual online: ned.is/hcbl400wt

Intended use

This product is intended to measure the diastolic and systolic blood pressure and pulse rate of an adult.

The product is not intended for use on neonates, children or pregnant women.

This product is not intended for medical use.

The product is intended for indoor use only.

Any modification of the product may have consequences for safety, warranty and proper functioning.

Specifications

Product	Blood Pressure Monitor
Article number	HCBL400WT
Dimensions (l x w x h)	136 x 120 x 61 mm
Weight (excl. batteries and cuff)	225g
Cuff circumference	22 - 42 cm
Measuring method	Oscillometric method, automatic inflation and measurement
Measurement range	
Cuff pressure	0 - 300 mmHg
Systolic	60 - 260 mmHg
Diastolic	40 - 199 mmHg
Pulse rate	40 - 180 beats/minute
Accuracy	
Pressure	± 3 mmHg
Pulse rate	Less than 60: ±3bpm More than 60 (incl.) : ±5%
Precision	1 mmHg
Memory	2 x 60 measurements
Battery type	4 x 1.5 V AA batteries (not included)
Battery life	Approximately 270 times
Operating temperature	10 - 40 °C
Storage temperature	-20 - 50 °C
Relative humidity	≤85 %
Environmental pressure	80 - 105 kPa

Main parts (image A)

- | | |
|---------------------------------|--|
| 1 Storage bag | 9 Low battery indicator |
| 2 Cuff | 10 Pressure in the cuff is instable or there is a lot of remnant air in the cuff |
| 3 Rubber air tube | 11 Blood pressure level classification indicator |
| 4 Air tube connector | 12 Time and date display |
| 5 Battery compartment | 13 Systolic blood pressure |
| 6 Memory button | 14 Diastolic blood pressure |
| 7 Start button | 15 Pulse rate |
| 8 Irregular heartbeat indicator | |

Safety instructions

⚠ WARNING

- Ensure you have fully read and understood the instructions in this document before you install or use the product. Keep the packaging and this document for future reference.
- Only use the product as described in this document.
- Do not use the product if a part is damaged or defective. Replace a damaged or defective product immediately.
- Do not drop the product and avoid bumping.
- Do not expose the product to water or moisture.
- Only use the provided power cable.
- Disconnect the product from the power source and other equipment if problems occur.
- Do not unplug the product by pulling on the cable. Always grasp the plug and pull.
- This product may only be serviced by a qualified technician for maintenance to reduce the risk of electric shock.
- Consult with a doctor before use if you have any medical concerns.
- Do not use this product in a moving vehicle. This may result in inaccurate measurements.

Inserting the batteries (image B)

1. Open the battery compartment **A5**.
2. Place the batteries.
3. Close **A5**.

Setting the time and date

1. Press and hold the start button **A7** and the memory button **A6** simultaneously for two seconds.

The time format starts to blink. The product is now in time format adjustment mode.

- ⓘ If no button is pressed for 60 seconds, the product exits time and date adjustment mode.
2. Press **A6** to change the time format.
 3. Press **A7** to confirm.

The product is now in year adjustment mode.

4. Press **A6** to adjust the number.
5. Press **A7** to confirm.

The product is now in month adjustment mode.

6. Repeat steps 4 – 5 to set the month, date, hour and minutes.

Connecting the cuff

1. Insert the rubber air tube **A3** into the air tube connector **A4**.
- ⓘ The connector must be completely inserted to prevent air leakage during blood pressure measurements.

Applying the cuff

1. Place the cuff **A2** around a bare upper arm.
- ⓘ **A2** must be placed 1-2 cm above the elbow joint.
2. Choose an arm to measure on. Always measure on the same arm.
- ⓘ If **A2** is applied on the left arm, **A3** must be on the inside of the arm, aligned with the middle finger.
 - ⓘ If **A2** is applied on the right arm, **A3** must be at the side of the elbow.
 - ⓘ **A2** should be tight but comfortable. It should be possible to insert a finger between the arm and **A2**.

Body posture during measurement

Sitting down

1. Sit down.
 2. Put your feet flat on the floor.
 3. Place your hand on a flat surface with the palm up.
- ⓘ Do not rest your arm on **A3**. This may restrict the airflow to **A2**.
4. Position your arm so **A2** is at the same height as your heart.
 5. Sit calmly for 5 minutes before starting the blood pressure measurement.

Lying down

1. Lie on your back.
2. Place your arm straight along your side with the palm up.
3. Position your arm so **A2** is at the same height as your heart.
4. Lie calmly for 5 minutes before starting the blood pressure measurement.

Measuring blood pressure

1. Press **A7**.
- ⓘ All symbols - **A8** to **A15** - are shown on the display. If not all symbols are visible, contact the service center.

The current memory group starts to blink in the place of **A14**.

- ⓘ Memory groups can be used to store measurements for multiple people on the same device.
2. Press **A6** to switch memory groups.
 3. Press **A7** to confirm.

- ⓘ The memory group is also confirmed automatically after 5 seconds of inactivity.

The display shows the most recent measurements.

If no measurements are stored in the memory group, the display shows zeroes.

The product starts to seek zero pressure.

A2 inflates until sufficient pressure has built up for a measurement.

A2 slowly deflates while the measurement is carried out.

The results are shown on the display.

- ⓘ The measurements are automatically stored.
4. Press **A7** to switch off the display manually.
- ⓘ The display switches off automatically after 1 minute of inactivity.

Reading measurements from the memory

1. Press **A6**.

The current memory group is shown in the place of **A14**.

2. Press **A7** to switch memory groups.
3. Press **A6** to confirm.

- ⓘ The memory group is also confirmed automatically after 5 seconds of inactivity.

The display shows the average values of all measurements in the current memory group.

If no measurements are stored in the memory group, the display shows zeroes.

4. Press **A6** to show the average values of the last three measurements in the current memory group.
5. Press **A6** to show the number of stored measurements in the current memory group.
6. Press **A6** to show the most recent measurements.
7. Press **A6** repeatedly to show the previous measurements.
8. Press **A7** to switch off the display manually.

- ⓘ The display switches off automatically after 1 minute of inactivity.

Deleting measurements from the memory

1. When the display shows a measurement, press and hold **A6** for three seconds to delete all measurements in the current memory bank.

Technical alarm

If the measured blood pressure is outside the rated range, the display shows 'Hi' or 'Lo' as technical alarm.

1. Check if the measurement was performed correctly.
2. If the measurement was performed correctly, consult your physician.

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